

EBP05854 - IMMUNIZING PEPTIDE

SKU: EBP05854



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_009213.1; NP_001365216.1; NP_001365220.1
Immunogen	Peptide with sequence C-PNGERDYNIAITV, from the C Terminus of the protein sequence according to NP_009213.1; NP_001365216.1; NP_001365220.1.
Product Comments	The peptide was used in the production of EB05854 - Goat Anti-RNF13 / RZF (isoform 1 and 3) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PNGERDYNIAITV
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)