

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05859 - IMMUNIZING PEPTIDE

SKU: EBP05859

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_078884.2; NP_001243481.1; NP_001339449.1

Peptide with sequence C-DIPSLEPEGTSAS, from the C Terminus of the protein sequence according to NP_078884.2; Immunogen NP_001243481.1; NP_001339449.1.

The peptide was used in the production of EB05859 - Goat Anti-NEIL1 / NEH1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-DIPSLEPEGTSAS Sequence

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$