

EBP05866 - IMMUNIZING PEPTIDE

SKU: EBP05866



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_733749.1; NP_001139674.1; NP_001139675.1; NP_001139676.1; NP_001139677.1; NP_001139678.1; NP_001139679.1
Immunogen	Peptide with sequence C-NSAQRRLGDSSDL, from the C Terminus of the protein sequence according to NP_733749.1; NP_001139674.1; NP_001139675.1; NP_001139676.1; NP_001139677.1; NP_001139678.1; NP_001139679.1.
Product Comments	The peptide was used in the production of EB05866 - Goat Anti-RASGRP4 (isoform 3) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NSAQRRLGDSSDL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)