

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05867 - IMMUNIZING PEPTIDE

SKU: EBP05867

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_067629.2; NP_001003827.1; NP_569073.1; NP_001003819.1 ID

Immunogen Peptide with sequence CPAPMTLCPPSS, from the C Terminus of the protein sequence according to NP_067629.2; NP_001003827.1; NP_569073.1; NP_001003819.1.

The peptide was used in the production of EB05867 - Goat Anti-TRIM34 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

CPAPMTLCPPSS

Sequence

Instructions Ambient Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade