

EBP05871 - IMMUNIZING PEPTIDE

SKU: EBP05871



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_067076.2; NP_742013.1
Immunogen	Peptide with sequence C-NAAPLTIRPPTDWE, from the C Terminus of the protein sequence according to NP_067076.2; NP_742013.1.
Product Comments	The peptide was used in the production of EB05871 - Goat Anti-RNF23 / TRIM39 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NAAPLTIRPPTDWE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)