

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05884 - IMMUNIZING PEPTIDE

SKU: EBP05884

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_060469.4; AAH12077.1 ID

 $\label{eq:lower_$

The peptide was used in the production of EB05884 - Goat Anti-RNF31 / ZIBRA (N terminus) Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

PGEEERAFLVAR-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade