

EBP05887 - IMMUNIZING PEPTIDE

SKU: EBP05887



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_055563

Immunogen Peptide with sequence HFSIPETESRSGD-C, from the N Terminus of the protein sequence according to NP_055563.

Product Comments The peptide was used in the production of EBP05887 - Goat Anti-SNX17 Antibody. It is currently untested in blocking applications.

Peptide Sequence HFSIPETESRSGD-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png