

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05904 - IMMUNIZING PEPTIDE**

**SKU:** EBP05904

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 071416.2; NP 690049.1; NP 690050.1 ID

Peptide with sequence ATPYVPVPMPIGN-C, from the N Terminus of the protein sequence according to Immunogen

NP 071416.2; NP 690049.1; NP 690050.1.

The peptide was used in the production of EB05904 - Goat Anti-SNX16 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** ATPYVPVPMPIGN-C Sequence

Shipping

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade