

# EBP05913 - IMMUNIZING PEPTIDE

**SKU:** EBP05913



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_150232

**Immunogen** Peptide with sequence PATPSLKVVHELPAC, from the N Terminus of the protein sequence according to NP\_150232.

**Product Comments** The peptide was used in the production of EBP05913 - Goat Anti-TRIM15 Antibody. It is currently untested in blocking applications.

**Peptide Sequence** PATPSLKVVHELPAC

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)

## DOCUMENTS

- [Data Sheet](#)