

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05914 - IMMUNIZING PEPTIDE

SKU: EBP05914

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Accession

Instructions Shipped at ambient temperature, store at -20°C.

ID

 $NP_056086.2;\ NP_001123539.1;\ NP_001337983.1;\ NP_001337984.1;\ NP_001337985.1l;\ NP_001337986.1$

Peptide with sequence DSGNHCFKVYRYLQ, from the C Terminus of the protein sequence according to NP_056086.2; Immunogen NP_001123539.1; NP_001337983.1; NP_001337984.1; NP_001337985.11; NP_001337986.1.

The peptide was used in the production of EB05914 - Goat Anti-TRIM2 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

DSGNHCFKVYRYLQ Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade