

EBP05929 - IMMUNIZING PEPTIDE

SKU: EBP05929



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_660095.1; NP_071763.2
Immunogen	Peptide with sequence GQESSKPVWPNPT-C, from the N Terminus of the protein sequence according to NP_660095.1; NP_071763.2.
Product Comments	The peptide was used in the production of EB05929 - Goat Anti-PJA1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	GQESSKPVWPNPT-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)