

EBP05943 - IMMUNIZING PEPTIDE

SKU: EBP05943



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001041666.1; NP_037414.3
Immunogen	Peptide with sequence C-LNQLFPGYGNGR, from the C Terminus of the protein sequence according to NP_001041666.1; NP_037414.3.
Product Comments	The peptide was used in the production of EB05943 - Goat Anti-ICBP90 / UHRF1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-LNQLFPGYGNGR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)