



EBP05944 - IMMUNIZING PEPTIDE

SKU: EBP05944



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_671725.1; NP_757379.1
Immunogen	Peptide with sequence C-QHAASLGLDTLPN, from the C Terminus of the protein sequence according to NP_671725.1; NP_757379.1.
Product Comments	The peptide was used in the production of EBP05944 - Goat Anti-DMBX1 / OTX3 / MBX Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QHAASLGLDTLPN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)