

EBP05945 - IMMUNIZING PEPTIDE

SKU: EBP05945



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003064.2; NP_001007469.1
Immunogen	Peptide with sequence MMMALSKTFGQKP-C, from the N Terminus of the protein sequence according to NP_003064.2; NP_001007469.1.
Product Comments	The peptide was used in the production of EBP05945 - Goat Anti-SMARCB1 / SNF5 / INI1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	MMMALSKTFGQKP-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiking_medium.png