

EBP05962 - IMMUNIZING PEPTIDE

SKU: EBP05962



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_116011

Immunogen Peptide with sequence C-RITLTIVSWPGTFL, from the C Terminus of the protein sequence according to NP_116011.

Product Comments The peptide was used in the production of EBP05962 - Goat Anti-LNX Antibody. It is currently untested in blocking applications.

Peptide Sequence C-RITLTIVSWPGTFL

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png