

EBP05971 - IMMUNIZING PEPTIDE

SKU: EBP05971



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055948.2; NP_001035529.1
Immunogen	Peptide with sequence MEIQMDEGGGVV-C, from the N Terminus of the protein sequence according to NP_055948.2; NP_001035529.1.
Product Comments	The peptide was used in the production of EB05971 - Goat Anti-JIP3 / Syd2 / JSAP1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	MEIQMDEGGGVV-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)