

EBP05986 - IMMUNIZING PEPTIDE

SKU: EBP05986



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_005142.1; NP_001032411.1

Immunogen Peptide with sequence PLYSVTVKWKGEK-C, from the N Terminus of the protein sequence according to NP_005142.1; NP_001032411.1.

Product Comments The peptide was used in the production of EBP05986 - Goat Anti-USP14 / TGT Antibody. It is currently untested in blocking applications.

Peptide Sequence PLYSVTVKWKGEK-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)