

# EBP05987 - IMMUNIZING PEPTIDE

**SKU:** EBP05987



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_115542.1, NP\_569055.1

**Immunogen** Peptide with sequence AGTSAPGSKRRSE-C, from the N Terminus of the protein sequence according to NP\_115542.1, NP\_569055.1.

**Product Comments** The peptide was used in the production of EBP05987 - Goat Anti-TREX1 / ATRIP (isoform) Antibody. It is currently untested in blocking applications.

**Peptide Sequence** AGTSAPGSKRRSE-C

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)