

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05995 - IMMUNIZING PEPTIDE

SKU: EBP05995

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_060568.3; NP_001245244.1; NP_001336023.1; NP_001336026.1; NP_001336030.1; NP_001336031.1; Accession NP_001336032.1 ID

Peptide with sequence AENSVLTSTTGRT-C, from the N Terminus of the protein sequence according to

Immunogen NP 060568.3; NP 001245244.1; NP 001336023.1; NP 001336026.1; NP 001336030.1; NP 001336031.1;

NP 001336032.1.

The peptide was used in the production of EB05995 - Goat Anti-ECT2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

AENSVLTSTTGRT-C Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade