



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05999 - IMMUNIZING PEPTIDE

SKU: EBP05999

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_004487.2 ID

Immunogen Peptide with sequence C-GVYSRPVLNTS, from the C Terminus of the protein sequence according to NP_004487.2.

Product The peptide was used in the production of EB05999 - Goat Anti-FOXA1 / HNF3A Antibody. It is currently untested in

Comments blocking applications.

Peptide

C-GVYSRPVLNTS Sequence

Shipping

Ambient Instructions

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$