

EBP06000-B - IMMUNIZING PEPTIDE

SKU: EBP06000-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_068556.2; NP_710141.1;
Immunogen	Peptide with sequence C-GVYSRPIMNSS., from the C Terminus of the protein sequence according to NP_068556.2; NP_710141.1;.
Product Comments	The peptide was used in the production of EB06000-B - Goat Anti-FOXA2 / HNF3B, Biotinylated Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GVYSRPIMNSS.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)