

EBP06026 - IMMUNIZING PEPTIDE

SKU: EBP06026



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060215.4; NP_001116855.1; NP_001116856.1; NP_001116857.1
Immunogen	Peptide with sequence LSATPLYGNVHSW-C, from the N Terminus of the protein sequence according to NP_060215.4; NP_001116855.1; NP_001116856.1; NP_001116857.1.
Product Comments	The peptide was used in the production of EBP06026 - Goat Anti-BCOR / BCL6 co-repressor Antibody. It is currently untested in blocking applications.
Peptide Sequence	LSATPLYGNVHSW-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)