

EBP06027 - IMMUNIZING PEPTIDE

SKU: EBP06027



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_001078940.1; NP_005188.2

Immunogen Peptide with sequence C-TAKGQKEQKETTKN, from the C Terminus of the protein sequence according to NP_001078940.1; NP_005188.2.

Product Comments The peptide was used in the production of EB06027 - Goat Anti-FOXN3 / CHES1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-TAKGQKEQKETTKN

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png