

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06030 - IMMUNIZING PEPTIDE

SKU: EBP06030

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001032208.1; NP_055729.2 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence C-NRNNSSGRVPTII, from the C Terminus of the protein sequence according to NP_001032208.1; NP_055729.2.}$

The peptide was used in the production of EB06030 - Goat Anti-CENTG2 / AGAP1 Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

C-NRNNSSGRVPTII Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade