

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06034 - IMMUNIZING PEPTIDE

SKU: EBP06034

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_004505.2; ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \ \ \ \ \ \ Peptide \ with \ sequence \ C-TPPAAVREKGVQN, from the C Terminus of the protein sequence according to NP_004505.2;.$

The peptide was used in the production of EB06034 - Goat Anti-FOXK2 / ILF (isoform 1) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-TPPAAVREKGVQN Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade