

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06046 - IMMUNIZING PEPTIDE

SKU: EBP06046

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_036294.2; NP_078831.4 ID

Immunogen Peptide with sequence CLEQLLTSPSPS, from the C Terminus of the protein sequence according to NP_036294.2; NP_078831.4.

The peptide was used in the production of EB06046 - Goat Anti-FBXL6 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

CLEQLLTSPSPS

Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade