

EBP06063 - IMMUNIZING PEPTIDE

SKU: EBP06063



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_009136; NP_542430; NP_542431; NP_542432
Immunogen	Peptide with sequence C-PMYLPDDPSLEA, from the C Terminus of the protein sequence according to NP_009136; NP_542430; NP_542431; NP_542432.
Product Comments	The peptide was used in the production of EBP06063 - Goat Anti-TREX2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PMYLPDDPSLEA
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png