

# EBP06079-B - IMMUNIZING PEPTIDE

**SKU:** EBP06079-B



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_006782.1; NP_996670.1
<b>Immunogen</b>	Peptide with sequence APKQDPKPKFQEGC., from the N Terminus of the protein sequence according to NP_006782.1; NP_996670.1.
<b>Product Comments</b>	The peptide was used in the production of EB06079-B - Goat Anti-MRG15 / MORF4L1, Biotinylated Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	APKQDPKPKFQEGC.
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>