

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06079 - IMMUNIZING PEPTIDE**

**SKU:** EBP06079

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 006782.1; NP 996670.1 ID

Peptide with sequence APKQDPKPKFQEG-C, from the N Terminus of the protein sequence according to **Immunogen** 

NP 006782.1; NP 996670.1.

The peptide was used in the production of EB06079 - Goat Anti-MRG15 / MORF4L1 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

APKQDPKPKFQEG-C Sequence

Shipping Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade