

EBP06079 - IMMUNIZING PEPTIDE

SKU: EBP06079



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006782.1; NP_996670.1
Immunogen	Peptide with sequence APKQDPKPKFQEG-C, from the N Terminus of the protein sequence according to NP_006782.1; NP_996670.1.
Product Comments	The peptide was used in the production of EB06079 - Goat Anti-MRG15 / MORF4L1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	APKQDPKPKFQEG-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)