

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06086 - IMMUNIZING PEPTIDE

SKU: EBP06086

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_056386.2; NP_001093879.1; NP_001291721.1 ID

Peptide with sequence AAGKSGGSAGEIT-C, from the N Terminus of the protein sequence according to

NP_056386.2; NP_001093879.1; NP_001291721.1.

The peptide was used in the production of EB06086 - Goat Anti-SUSP1 / SSP1 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

AAGKSGGSAGEIT-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade