

EBP06088 - IMMUNIZING PEPTIDE

SKU: EBP06088



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_002810.1; NP_114367.1; NP_114368.1; NP_787041.1
Immunogen	Peptide with sequence DGIVPDIAVGTKR-C, from the N Terminus of the protein sequence according to NP_002810.1; NP_114367.1; NP_114368.1; NP_787041.1.
Product Comments	The peptide was used in the production of EBP06088 - Goat Anti-PTBP1 / PTB Antibody. It is currently untested in blocking applications.
Peptide Sequence	DGIVPDIAVGTKR-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)