

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06088 - IMMUNIZING PEPTIDE

SKU: EBP06088

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

ID

NP_002810.1; NP_114367.1; NP_114368.1; NP_787041.1

Peptide with sequence DGIVPDIAVGTKR-C, from the N Terminus of the protein sequence according to Immunogen

NP 002810.1; NP 114367.1; NP 114368.1; NP 787041.1.

The peptide was used in the production of EB06088 - Goat Anti-PTBP1 / PTB Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

DGIVPDIAVGTKR-C Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png