

# EBP06090 - IMMUNIZING PEPTIDE

**SKU:** EBP06090



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_065705

**Immunogen** Peptide with sequence DKRKLGRPPSSSE-C, from the N Terminus of the protein sequence according to NP\_065705.

**Product Comments** The peptide was used in the production of EBP06090 - Goat Anti-SEN7 (N Terminus) Antibody. It is currently untested in blocking applications.

**Peptide Sequence** DKRKLGRPPSSSE-C

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)