



EBP06099 - IMMUNIZING PEPTIDE

SKU: EBP06099



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_055647.2

Immunogen Peptide with sequence CNPNNKAKIGNKP, from the C Terminus of the protein sequence according to NP_055647.2.

Product Comments The peptide was used in the production of EBP06099 - Goat Anti-TBC1D4 Antibody. It is currently untested in blocking applications.

Peptide Sequence CNPNNKAKIGNKP

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png