

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06119 - IMMUNIZING PEPTIDE**

**SKU:** EBP06119

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_008920.5; NP\_001073346.2 ID

Immunogen Peptide with sequence AKRSSLYIRIVEG-C, from the N Terminus of the protein sequence according to NP\_008920.5; NP\_001073346.2.

The peptide was used in the production of EB06119 - Goat Anti-CAPRI / RASA4 Antibody. It is currently untested in **Product** 

**Comments** blocking applications.

**Peptide** 

AKRSSLYIRIVEG-C Sequence

Ambient Instructions Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade