

EBP06119 - IMMUNIZING PEPTIDE

SKU: EBP06119



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_008920.5; NP_001073346.2
Immunogen	Peptide with sequence AKRSSLYIRIVEG-C, from the N Terminus of the protein sequence according to NP_008920.5; NP_001073346.2.
Product Comments	The peptide was used in the production of EB06119 - Goat Anti-CAPRI / RASA4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	AKRSSLYIRIVEG-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)