

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06126 - IMMUNIZING PEPTIDE

SKU: EBP06126

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 056057.2; NP 001035207.1; NP 001128662.1 ID

Peptide with sequence C-AFTADPLSLLRNV, from the C Terminus of the protein sequence according to Immunogen

 $NP_056057.2;\ NP_001035207.1;\ NP_001128662.1.$

The peptide was used in the production of EB06126 - Goat Anti-CENTD2 / ARAP1 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-AFTADPLSLLRNV Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade