

EBP06133 - IMMUNIZING PEPTIDE

SKU: EBP06133



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_056125.3; NP_996792.1
Immunogen	Peptide with sequence C-DHRPETPSIKQKVT, from the C Terminus of the protein sequence according to NP_056125.3; NP_996792.1.
Product Comments	The peptide was used in the production of EB06133 - Goat Anti-HCA67 / PSD3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DHRPETPSIKQKVT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)