

EBP06138 - IMMUNIZING PEPTIDE

SKU: EBP06138



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_054865.2; NP_115852.1; NP_115853.2
Immunogen	Peptide with sequence C-PDYVNGEVAATEA, from the C Terminus of the protein sequence according to NP_054865.2; NP_115852.1; NP_115853.2.
Product Comments	The peptide was used in the production of EB06138 - Goat Anti-LAT2 / NTAL / LAB Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PDYVNGEVAATEA
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png