

EBP06142 - IMMUNIZING PEPTIDE

SKU: EBP06142



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_055886.1; NP_001129080.1
Immunogen	Peptide with sequence C-KTGLIPENYVEFL, from the C Terminus of the protein sequence according to NP_055886.1; NP_001129080.1.
Product	The peptide was used in the production of EB06142 - Goat Anti-GRAF / ARHGAP26 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-KTGLIPENYVEFL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)