

EBP06150 - IMMUNIZING PEPTIDE

SKU: EBP06150



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_110385.1; NP_001269667.1; NP_001269668.1; NP_001269669.1; NP_001269670.1; NP_001269671.1
Immunogen	Peptide with sequence KEKQDWHMCPNIF, from the C Terminus of the protein sequence according to NP_110385.1; NP_001269667.1; NP_001269668.1; NP_001269669.1; NP_001269670.1; NP_001269671.1.
Product Comments	The peptide was used in the production of EB06150 - Goat Anti-OSBP2 / ORP4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	KEKQDWHMCPNIF
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)