

EBP06153 - IMMUNIZING PEPTIDE

SKU: EBP06153

SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_056365.1; NP_663160.1; NP_663161.1; NP_663162.1
Immunogen	Peptide with sequence MSDEKNLGVSQKL-C, from the N Terminus of the protein sequence according to NP_056365.1; NP_663160.1; NP_663161.1; NP_663162.1.
Product Comments	The peptide was used in the production of EB06153 - Goat Anti-ORP3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	MSDEKNLGVSQKL-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)