

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06154 - IMMUNIZING PEPTIDE

SKU: EBP06154

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_065947.1; NP_663613.1; NP_001137535.1 ID

Peptide with sequence KEEAFLRRRFSLC, from the N Terminus of the protein sequence according to NP_065947.1; Immunogen NP_663613.1; NP_001137535.1.

The peptide was used in the production of EB06154 - Goat Anti-ORP5 (OBPH1) Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

KEEAFLRRRFSLC Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade