

EBP06155 - IMMUNIZING PEPTIDE

SKU: EBP06155

SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_115912.1; NP_665682.1; NP_001188409.1; NP_001188410.1; NP_001188411.1
Immunogen	Peptide with sequence C-DPGFSKVDSPVLW, from the C Terminus of the protein sequence according to NP_115912.1; NP_665682.1; NP_001188409.1; NP_001188410.1; NP_001188411.1.
Product Comments	The peptide was used in the production of EBP06155 - Goat Anti-OSBPL6 / ORP6 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DPGFSKVDSPVLW
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)