

EBP06158 - IMMUNIZING PEPTIDE

SKU: EBP06158



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_078862.4; NP_683702.1; NP_683704.2; NP_683705.1; NP_683706.3; NP_683707.3
Immunogen	Peptide with sequence C-EPLLKRLGAAKH, from the C Terminus of the protein sequence according to NP_078862.4; NP_683702.1; NP_683704.2; NP_683705.1; NP_683706.3; NP_683707.3.
Product Comments	The peptide was used in the production of EB06158 - Goat Anti-ORP9 / OSBPL9 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EPLLKRLGAAKH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)