

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06158 - IMMUNIZING PEPTIDE**

**SKU:** EBP06158

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 078862.4; NP 683702.1; NP 683704.2; NP 683705.1; NP 683706.3; NP 683707.3 ID

Peptide with sequence C-EPLLKRLGAAKH, from the C Terminus of the protein sequence according to Immunogen

 $NP\_078862.4; \ NP\_683702.1; \ NP\_683704.2; \ NP\_683705.1; \ NP\_683706.3; \ NP\_683707.3.$ 

The peptide was used in the production of EB06158 - Goat Anti-ORP9 / OSBPL9 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-EPLLKRLGAAKH Sequence

**Shipping** 

Grade

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png