

EBP06161 - IMMUNIZING PEPTIDE

SKU: EBP06161



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004811.1
Immunogen	Peptide with sequence C-YPDSDVLFYKVKKS, from the C Terminus of the protein sequence according to NP_004811.1.
Product Comments	The peptide was used in the production of EB06161 - Goat Anti-CYP7B1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-YPDSDVLFYKVKKS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)