

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06170 - IMMUNIZING PEPTIDE

SKU: EBP06170

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_057529 ID

Immunogen Peptide with sequence C-ERAEWIEAIKKLT, from the C Terminus of the protein sequence according to NP_057529.

Product The peptide was used in the production of EB06170 - Goat Anti-Pleckstrin 2 Antibody. It is currently untested in

Comments blocking applications.

Peptide

Sequence

C-ERAEWIEAIKKLT

Shipping Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png