

EBP06171 - IMMUNIZING PEPTIDE

SKU: EBP06171



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_037517.1

Immunogen Peptide with sequence DLCHPEPAELSSG, from the N Terminus of the protein sequence according to NP_037517.1.

Product Comments The peptide was used in the production of EB06171 - Goat Anti-PSCD4 / CYTOHESIN 4 Antibody. It is currently untested in blocking applications.

Peptide Sequence DLCHPEPAELSSG

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)